

Table. Assays Used at CHECKMATE Sites and Their Diagnostic Cut-points

Site	CK-MB Assay	CK-MB (mcg/L)	Troponin Assay	TnI (mcg/L)	TnT (mcg/L)
1	Abbott Axsym	10.0	Abbott Axsym	1.0	---
2	Roche Elecsys	10.0	Roche Elecsys	---	0.1
3	Abbott Axsym	9.0	Roche Elecsys	---	0.1
4	Abbott Axsym	4.0	---	---	---
5	Dade-Behring Stratus CS	5.1	Dade-Behring Stratus CS	0.4	---
6	Dade-Behring Opus Magnum	5.0	Dade-Behring Opus Magnum	0.5	---